

Musician ER Earmolds




Oto Seal



- Beige (default)
- Available in all colors
- Countersunk if requested

Silicone



- Opaque in color (default)
- Available in all colors
- Countersunk only

Musician ER Earmolds


High-Fidelity Hearing Protection

Music • Clubs • Concerts

Preserve the timbre and richness of music!

Musician Earplugs were the first high fidelity hearing protection available to musicians. OTO Hearing Products is an authorized manufacturer by Etymotic Research to produce ER Musician ear plugs. We manufacture these in two different materials; OTO Seal which is a rigid material that will become semi-flexible from body temperature when worn or a soft high grade silicone (illustrated above).

Musicians practice and perform in a variety of different settings and are exposed to high levels of sound, sometimes for long periods of time. They require different amounts of protection, depending on the sound levels they encounter during rehearsals and performances. Musician Earplugs replicate the ear canals natural response so that the sound heard with these earplugs has the same quality as the original, just more quiet. As designated by the number classification, ER 25's and ER 9's provide either more or less attenuation than ER 15's. Some musicians use one type of attenuator in one ear and one type in the other ear, depending on the source and location of the sound.



Three different types of interchangeable filters provide various types of sound reduction.

www.otohearing.com

Oto Hearing Products
1-800-661-7723 Toll Free
604-590-1211 Local